

Endo-Boot™

0071**1621**

3-7mm inner diameter,
4cm outer diameter
14cm length
for use on all scopes and delicate instruments
up to 5mm in diameter

0071**1622**

7-18mm inner diameter,
4cm outer diameter
14cm length
for use on all scopes and delicate
instruments over 5mm in diameter



The Endo-Boot™ helps protect delicate endoscopes from accidental stretching and tearing of the bending rubber section, as well as chipping or cracking of the distal lens. Facilities that use the Endo-Boot™ have reported a 50% reduction in bending rubber section replacements¹ – vastly reducing repair costs and downtime of instruments.

1621



1622

simple, easy-to-use design

*inexpensive, smart insurance cuts costly
and timely repairs*

smart protection

Source¹: Gastroenterology Nursing 11/2000
Endo-Boot™ is a trademark of US Endoscopy