


single-use

self-expanding Duodenal & Colonic Stents



Characteristics

- TTS (Through-The-Scope)
- Application Set Fr. 10
- Length:
 - 210 cm (duodenal)**
 - 230 cm (colonic)**
- Recommended guide wire: .035 inch (0.9 mm)

Placement method

- Distal


Special features

- braided Teflon tubing with transparent distal end
- easy handling and optimal deployment time
- no stent shortening at distal end of release
- stent design with smooth atraumatic stent ends
- effective anchorage and migration prevention
- constant and atraumatic structure of the braided design at any angle



self-expanding Duodenal Stents

- incl. pre-loaded application set
- distal placement
- radiopaque markers




uncoated

Order No.	length	stent diameter	tulip diameter
DUO-0-20-60	60 mm	20 mm	24 mm
DUO-0-20-80	80 mm	20 mm	24 mm
DUO-0-20-100	100 mm	20 mm	24 mm

self-expanding Colonic Stents

- incl. pre-loaded application set
- distal placement
- radiopaque markers



uncoated

Order No.	length	stent diameter	tulip diameter
COL-0-20-80	80 mm	20 mm	24 mm
COL-0-20-100	100 mm	20 mm	24 mm
COL-0-20-120	120 mm	20 mm	24 mm
COL-0-22-80	80 mm	22 mm	26 mm
COL-0-22-100	100 mm	22 mm	26 mm
COL-0-22-120	120 mm	22 mm	26 mm