

Basket Sizes

ENDO-FLEX baskets sizes are measured in stretched condition, which means when the basket is retracted into the tube. Other manufacturers may measure basket sizes in opened condition and by indicating the height and width in mm.

Conversion table for the different measuring methods:


ENDO-FLEX (stretched condition)	others (height x width)
40 mm	~ 28 x 15
50 mm	~ 37 x 20
60 mm	~ 45 x 25
70 mm	~ 50 x 30

Ø 2.3 mm

134 °C
273 F

4 wires

- detachable for cleaning



dormia-type

basket size (in stretched condition)

40 mm 50 mm 60 mm 70 mm

basket-type


basket size (in stretched condition)

40 mm 50 mm 60 mm 70 mm

	40 mm	50 mm	60 mm	70 mm	40 mm	50 mm	60 mm	70 mm
monofilament wire	150 224	150 225	150 226	150 227	151 224	151 225	151 226	151 227
braided wire	160 224	160 225	160 226	160 227	161 224	161 225	161 226	161 227

6 wires

- detachable for cleaning



dormia-type

basket size (in stretched condition)

40 mm 50 mm 60 mm 70 mm

basket-type


basket size (in stretched condition)

40 mm 50 mm 60 mm 70 mm

	40 mm	50 mm	60 mm	70 mm	40 mm	50 mm	60 mm	70 mm
monofilament wire	152 224	152 225	152 226	152 227	153 224	153 225	153 226	153 227
braided wire	162 224	162 225	162 226	162 227	163 224	163 225	163 226	163 227

4 wires / Guide wire

- accepts guide wire .035 inch



dormia-type

basket size (in stretched condition)

40 mm 50 mm 60 mm 70 mm

basket-type

basket size (in stretched condition)

40 mm 50 mm 60 mm 70 mm

	40 mm	50 mm	60 mm	70 mm	40 mm	50 mm	60 mm	70 mm
monofilament wire	150 224GW	150 225GW	150 226GW	150 227GW	151 224GW	151 225GW	151 226GW	151 227GW
braided wire	160 224GW	160 225GW	160 226GW	160 227GW	161 224GW	161 225GW	161 226GW	161 227GW




Ø 2.3 mm

134 °C
273 °F

6 wires / Guide wire

- accepts guide wire .035 inch



dormia-type

basket size (in stretched condition)

40 mm 50 mm 60 mm 70 mm

basket-type

basket size (in stretched condition)


40 mm 50 mm 60 mm 70 mm

monofilament wire	152 224GW	152 225GW	152 226GW	152 227GW	153 224GW	153 225GW	153 226GW	153 227GW
-------------------	-----------	-----------	-----------	-----------	-----------	-----------	-----------	-----------

Ø 1.8 mm

4 wires

- detachable for cleaning



dormia-type

basket size (in stretched condition)

40 mm 50 mm 60 mm 70 mm

basket-type


basket size (in stretched condition)

40 mm 50 mm 60 mm 70 mm

monofilament wire	150 184	150 185	150 186	150 187	151 184	151 185	151 186	151 187
braided wire	160 184	160 185	160 186	160 187	161 184	161 185	161 186	161 187

6 wires

- detachable for cleaning



dormia-type

basket size (in stretched condition)

40 mm 50 mm 60 mm 70 mm

basket-type

basket size (in stretched condition)

40 mm 50 mm 60 mm 70 mm

monofilament wire	152 184	152 185	152 186	152 187	153 184	153 185	153 186	153 187
-------------------	---------	---------	---------	---------	---------	---------	---------	---------



Ø 2.3 mm

134 °C
273 F

4 wires, Nitinol

- basket made of gold plated Nitinol, radiopaque tube, incl. handle



dormia-type

basket size (in stretched condition)

40 mm 50 mm 60 mm

basket-type

basket size (in stretched condition)

40 mm 50 mm 60 mm

Nitinol wire	146 224	146 225	146 226	141 224	141 225	141 226
--------------	---------	---------	---------	---------	---------	---------

6 wires, Nitinol

- basket made of gold plated Nitinol, radiopaque tube, incl. handle



dormia-type

basket size (in stretched condition)

40 mm 50 mm 60 mm

basket-type

basket size (in stretched condition)

40 mm 50 mm 60 mm



Nitinol wire	147 224	147 225	147 226	142 224	142 225	142 226
--------------	---------	---------	---------	---------	---------	---------

Ø 2.6 mm

134 °C
273 F

4 wires, Nitinol

- basket made of gold plated Nitinol, radiopaque tube, incl. handle
- guide wire compatible .035 inch



dormia-type

basket size (in stretched condition)

40 mm 50 mm 60 mm

basket-type

basket size (in stretched condition)

40 mm 50 mm 60 mm

monofilament wire	146 264GW	146 265GW	146 266GW	141 264GW	141 265GW	141 266GW
-------------------	-----------	-----------	-----------	-----------	-----------	-----------



Ø 1.8 mm

134 °C
273 °F

4 wires, Nitinol

- basket made of gold plated Nitinol, radiopaque tube, incl. handle



dormia-type

basket size (in stretched condition)

40 mm 50 mm 60 mm



basket-type

basket size (in stretched condition)

40 mm 50 mm 60 mm

Nitinol wire	146 184	146 185	146 186	141 184	141 185	141 186
--------------	---------	---------	---------	---------	---------	---------

6 wires, Nitinol

- basket made of gold plated Nitinol, radiopaque tube, incl. handle



dormia-type

basket size (in stretched condition)

40 mm 50 mm 60 mm



basket-type

basket size (in stretched condition)

40 mm 50 mm 60 mm

Nitinol wire	147 184	147 185	147 186	142 184	142 185	142 186
--------------	---------	---------	---------	---------	---------	---------

