

The Dream Team...

of foreign body management devices



The Dream Team...

US Endoscopy has redefined the standard of care in foreign body management by creating innovative and practical devices that specifically address these challenging procedures.



tapered tip

Guardus® overtube

Gentle protection and enhanced safety during foreign body extraction.

BX00711146 / BXSI711146 (25cm length)
BX00711148 / BXSI711148 (50cm length)
standard scopes (8.6-10.0mm scope OD)

BX00711147 / BXSI711147 (25cm length)
BX00711149 / BXSI711149 (50cm length)
therapeutic scopes (10.0-11.7mm scope OD)

enhanced
safety

Raptor® grasping device

Offers the versatility needed for retrieving a variety of foreign bodies and stents.

BX00711177
sheath diameter 2.4mm
length 230cm

NEW
Raptor® grasping device-mini
BX00711178
sheath diameter 1.8mm
length 230cm



Talon® grasping device

Grips into the food impaction for a swift retrieval.

BX00711175
sheath diameter 2.5mm
length 160cm

Roth Net® retrievers

From small to large foreign bodies to food bolus removal, the Roth Net® retrievers are designed to fully encapsulate objects for secure retrieval.

Roth Net® retriever - standard

BX00711050 /
BX00711162*
net size 3.0cm x 6.0cm
sheath diameter 2.5mm
length 230cm



Roth Net® retriever - mini

BX00711057
net size 2.0cm x 4.5cm
sheath diameter 1.8mm
length 160cm

NEW
BX00711168*
net size 2.0cm x 4.5cm
sheath diameter 1.8mm
length 160cm



Roth Net® Platinum® retriever - universal

BX00715050 /
BX00715060*
net size 4.0cm x 5.5cm
sheath diameter 2.5mm
length 230cm



Roth Net® Platinum® retriever - food bolus

BX00711155
net size 4.0cm x 5.5cm
sheath diameter 2.5mm
length 160cm



Roth Net® retriever - maxi

BX00711059
net size 4.0cm x 8.0cm
sheath diameter 3.0mm
length 160cm

BX00711165*
net size 4.0cm x 8.0cm
sheath diameter 3.0mm
length 230cm



Falcon® rotatable retrieval basket

A great utility tool to have on hand for removal of long objects and other foreign bodies not suitable for the Roth Net® retriever.

BX00711141
basket size 3.5cm x 6.0cm
sheath diameter 2.5mm
length 230cm



*sterile



listening... and delivering solutions®

5976 Heisley Road
Mentor, OH 44060
phone (1) 440 / 639-4494
fax (1) 440 / 639-4495
customer service (1) 800 / 769-8226
email global@usendoscopy.com

www.usendoscopy.com