



STERILE EO

Biliary & Esophageal / Colonic Multi-Stage Dilation Balloons

- for guide wire .035 inch
- radiopaque markers
- 3 progressive stages



Biliary Dilation Balloons	balloon-Ø (under pressure of)			balloon length	catheter length	catheter-Ø
3410620PRO	6 mm (3.0 atm)	7 mm (6,0 atm)	8 mm (9.0 atm)	20 mm	200 cm	7 Fr.
3410830PRO	8 mm (3.0 atm)	9 mm (5.5 atm)	10 mm (8.0 atm)	30 mm	200 cm	7 Fr.
3411030PRO	10 mm (3.0 atm)	11 mm (5.0 atm)	12 mm (7.0 atm)	30 mm	200 cm	7 Fr.
3411230PRO	12 mm (3.0 atm)	13.5 mm (4.5 atm)	15 mm (7.5 atm)	30 mm	200 cm	7 Fr.

Esophageal / Colonic Dilation Balloons	balloon-Ø (under pressure of)			balloon length	catheter length	catheter-Ø
34106PRO	6 mm (3.0 atm)	7 mm (6,0 atm)	8 mm (9.0 atm)	55 mm	230 cm	7 Fr.
34108PRO	8 mm (3.0 atm)	9 mm (5.5 atm)	10 mm (8.0 atm)	55 mm	230 cm	7 Fr.
34110PRO	10 mm (3.0 atm)	11 mm (5.0 atm)	12 mm (7.0 atm)	55 mm	230 cm	7 Fr.
34112PRO	12 mm (3.0 atm)	13.5 mm (4.5 atm)	15 mm (7.5 atm)	55 mm	230 cm	7 Fr.
34115PRO	15 mm (3.0 atm)	16.5 mm (4.5 atm)	18 mm (6.5 atm)	55 mm	230 cm	7 Fr.
34118PRO	18 mm (3.0 atm)	19 mm (4.5 atm)	20 mm (5.5 atm)	55 mm	230 cm	7 Fr.

Two types of applications combined in one balloon:

1. Direct application by leaving the inserted stylet / guide wire in place.
2. OTW application over any guide wire .035" after removing the stylet.

recommended equipment

Inflation / Deflation Device
incl. PSI-Manometer 0-440 PSI



Order No.

340 800 INFL

