




STERILE EO

Single-Stage-Balloon



triple lumen



air capacity (ml)	Ø balloon inflated (mm)	catheter	fitting to guide wire (inch)	Order No.
5.0	16.0	Fr. 7.0	.035	134 030
4.0	14.0	Fr. 7.0	.035	134 031
2.4	12.0	Fr. 7.0	.035	134 032
1.7	9.0	Fr. 7.0	.035	134 033

double lumen



air capacity (ml)	Ø balloon inflated (mm)	catheter	fitting to guide wire (inch)	Order No.
1.5	9.0	Fr. 5.0	.025	134 000
2.5	12.0	Fr. 5.0	.025	134 005
2.5	14.0	Fr. 7.0	.035	134 010
2.5	14.0	Fr. 5.0 - Fr. 7.0	.035	134 011
4.0	18.0	Fr. 7.0	.035	134 020